Powered by Microsoft® Windows®CE



# MONITOR

Laser Power / Energy Meter

((

- Single Channel
- Microsoft Windows CE™
- User-Friendly
- Large Display
- Full Statistical Functions
- Power/Energy/Photo Detectors
- USB / RS-232

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and/or other countries.

#### SOLO

Gentec Electro-Optics SOLO is the newest, most advanced, microprocessor-based power and energy meter in the Gentec-EO arsenal. It is the ideal choice for applications that only need to read one detector at a time. We have added more capability and at the same time made our friendly interface even more user friendly. The SOLO gives you complete statistics and allows you to save data internally and to your PC.

#### Windows CE™

SOLO's menu-driven Windows CE™ based interface is easy and intuitive enough to master in minutes. You can also download the operating software for the SOLO from our website. Our newest features, latest improvements, and even custom modifications just for you, are available at the click of a button.

### Power and Energy

The SOLO PE gives you the flexibility to read both power and energy detectors as simply as plugging in a new head. The SOLO lets you display power from a photodiode detector in dBm or watts. The SOLO will recognize the detector and automatically download all the calibration information it needs from the detector's EEPROM.

#### **OEM Version**

SOLO X is our PC-based version of the SOLO PE, without keyboard or display. It is ideal for networked OEM installations. It is great where you have a high degree of automation or need to keep people away from the process.

# Large Display

Your eyes will appreciate the easy to view high resolution 58x38 mm graphic LCD display of the SOLO PE. Your computer screen is the display for the SOLO X.

# **Ergonomic Design**


Compact, sturdy, and easy to handle, the SOLO can go anywhere! You will appreciate its ergonomic shape and adjustable stand. The stand protects the SOLO when you travel and it is adjustable so you can avoid the glare from overhead lights.

#### Data Acquisition

Data acquisition and remote control are possible via serial RS-232 and USB ports for the SOLO PE, and via USB port for the SOLO X. Your can use an adapter to plug the SOLO into your GPIB system and have the same functionality. You can save your settings so the SOLO PE will be ready to go as soon as you power it up. You can turn on the time stamp to ease your analysis or leave it off to maximize your storage. Data is easily plotted on widely available spreadsheets and plotting programs.

#### Remote Control

PC-SOLO, a user-friendly communication software package, lets you control the SOLO PE and X from your PC. It is particularly useful for data acquisition since the storage capacity is determined by your PC's hard disk. It is the only interface for the SOLO X. Upgrades for PC-SOLO are available on our website. You'll find our LabVIEW drivers there as well. Use them as is, or customize them for your specific needs.




PC-SOLO Software



- 0 to 1 volt analog output
- Detector
- RS-232 connector
- External power supply input
- Mini USB connector

#### CONVENIENT WINDOWS CE™ MENUS



SOLO's menu-driven Windows CE™ based interface is easy and intuitive enough to master in minutes.

#### Settings Display Scale Acquire Data Analog Out Data Transfer ✓ Anticipation Statistics Mode > Display Hi Res Zero Offset Lock Keyboard Calo Mode ad=2.00mW Res=14.55uW @ 1064 nm


A multitude of easy window based controls for every application.

#### **REAL TIME DISPLAY**

#### Scale Settings Ctrl Display 15.20 mW 0.00 mW 30.00 mW Head=15.20mW Res=14.65uW @ 1064 nm


High or Low Resolution. Select the high resolution mode for the most significant digits available or the standard resolution to filter out unimportant fluctuations in the measurement.

### HISTOGRAM DISPLAY




Easily switch between Real time, Line plot, Histogram, Digital needle and Statistics displays.

#### STATISTICS DISPLAY



The SOLO can display a complete statistical analysis of power or energy measurements.

#### DATA SAMPLING SETUP




You have complete control over the data sampling. Use the defaults or select your own sample rate, sample period, and the time period or number of points to do the statistics.

#### TUNING DISPLAY



The SOLO also provides a fast, high resolution digital needle suitable for tuning and aligning. Refresh rate: 15Hz.

#### STATUS SCREEN



Look at the Status screen to review all of your settings at any time.

### **Automatic Head Recognition**

The SOLO works with all Gentec-EO detectors. It automatically recognizes every detector head, ensuring accurate auto-calibration. More importantly, it can take advantage of our Personal wavelength correction™. It reads the memory in the Smart Interface connector to provide a wavelength correction that is based on spectral data measured from that specific detector. Your measurements across the band have never been this precise or easy.

### Fast Response

The SOLO can handle data from energy detectors up to 1 kHz. It also features an anticipation circuit for fast power detector response.

## **Quality Thresholds**

We even provide threshold lines you can overlay on your line plot so you can easily tell when the measurement is within or outside of tolerance. This is excellent for quality control and pass/fail monitoring applications.

#### **Correction Factors**

If you use a sampling device such as the Gentec-EO Holographic Beam Sampler (HBS) or need to account for other optics in your system you can easily program correction factors to obtain a true value reading. The SOLO allows you to enter two multipliers and 2 offsets for maximum measurement flexibility.

### Saves Data and Settings

The SOLO can store data in flash memory and transfer it to a PC. The SOLO can also remember and recall your settings. Look at the Status screen to review all of your settings at any time.

#### Flexible Power

Use the battery in the field. Recharge and take data at the same time with the external power supply. In addition, operate the SOLO without battery or power supply when plugged into a USB port! This is great for changing batteries without interrupting your measurements.

#### Flexible Statistics

You can set the SOLO to calculate the statistics for a single sample and stop, or to repeat continuously. Take data for a few seconds or a few weeks. You have the flexibility to handle any application from analyzing a single short pulse with high resolution to sampling performance over a period of months.







PROTECTIVE

POUCH



SOLO X



OPTIONAL WALL BRACKET



TNCLLIDED LISB CARLE



OPTIONAL RS-232 CABLE



BATTERY PACK

	SOLO PE	SOLO X	
POWER METER SPECIFICATIONS			
Power Range (thermal head)	100 μW to 10 kW	✓	
Power Range (photo detector head)	1 pW to 3 W	<b>√</b>	
Resolution (digital)	15 μW on the 30 mW scale	✓	
Monitor Accuracy	± 1.0% a	✓	
hermopile Response Time (accelerated)			
Statistics	Current value, Max, Min, Average,		
	Std Dev., RMS stability, PTP stability, Time	<b>✓</b>	
NERGY METER SPECIFICATIONS			
inergy Range	8 μJ to 20 kJ	✓	
esolution (digital)	50 nJ	✓	
Ionitor Accuracy	1.0% a	✓	
oftware Trigger Level	0.1 to 99.9%, 0.1% resolution, default 2% <sup>a</sup>	✓	
Repetition Rate	1000 Hz	✓	
Statistics	Current value, Max, Min, Average, Std Dev., RMS stability, PTP stability, Pulse #, Repetition Rate, Avg Power	✓	
DETECTOR COMPATIBILITY			
Thermopile (power detector)	Average Power	✓	
,	Single Shot Energy ✓		
Pyroelectric (energy detector)	Pulse Energy including single shot	✓	
Photo Detector	Average Power (mW, dBm)	✓	
GENERAL SPECIFICATIONS			
Digital Display	58 x 38 mm LCD, 160 x 240 pixels	PC screen	
Data Storage	up to 225 000 points <sup>b</sup>	✓	
Data Displays	Real time, Line plot, Histogram,		
	Statistics, Tuning display	✓	
High Resolution mode	More Significant Digits	✓	
User input correction factors	2 Multipliers and 2 Offsets (7 digits)	✓ .	
Analog Output	0-1 volt, full scale, ±1%	✓	
Internet upgrades	USB and RS-232° (up to 115.2 kbps)	USB	
PC SOLO	✓	✓	
PC Serial Commands	USB and RS-232°	USB	
LabVIEW Drivers	✓	✓	
Dimensions (without stand)	230 (L) x 95 max (W) x 55 max (H) mm	231 x 102 x 35 mm	
Weight	0.58 kg (with stand)	0.44 kg	
Battery pack	4 rechargeable 1.2 V Ni-MH AA		
	(6 hours to full charge)		
External Power Supply	220/240 VAC, 50 Hz to 12 VDC, 1A;		
	100/120 VAC, 60 Hz to 12 VDC, 800 mA; or USB	USB	

Order number	Description	Order number	Description
100-200076	SOLO PE + Battery + Power Supply + USB Cable	100-200128	Optional SOLO Protective Pouch
100-200293	SOLO X + USB Cable	100-100283	Optional SOLO Wall Bracket
100-200106	Optional SOLO RS-232 Cable	100-200083 Spare SOLO Battery Pack	

a. For Rep. Rate < 500 Hz.

Specifications subject to change without notice.

b. Maximum when not using advanced features that require memory.

c. USB Cable is included with the SOLO, RS 232 cable must be purchased separately.





UP19K-15S-H5

UP19K-30H-H5

# POWER DETECTORS

Mid Power - Compact

 $\epsilon$ 

- Fast
- Small Size
- Flat Spectral Response
- Full NIST-Traceability
- High Damage Threshold: 36 kW/cm²
- Personal Wavelength Correction™

#### Ultra Series UP19K-H5

Ultra Power and Ultra Performance are what you get with our new UP series detectors. Ultra performance means fast. These are the fastest detectors in the Gentec-EO arsenal. Ultra performance means compact. The UP series detectors can take a lot of laser power in a small body. Ultra performance means flexible. They come ready to mount on a rod, a bracket and the square case even lets you set them right on the table. Ultra performance means expandable. We can easily increase the power capability of your modular UP series detector as your needs change. Ultra performance means accurate. It is hard to do better than our NIST traceable calibration and Personal wavelength correction™. Ultra performance means versatile. Each model is compatible with all Gentec-EO monitors. A UP series detector is the best choice for many applications.

#### Low Profile

The UP19K-15S-H5, just 21 mm thin, is ideal for laser maintenance and service applications. From a few mW to 15 W it is the smallest, toughest and fastest detector in its class. The absorber handles 36 kW/cm² of average power density. Only Gentec-E0's WB detectors can do better. It is fast with a response better than 0.6 seconds using a Gentec-E0 monitor. You can even run it up to 23 watts for short periods. That's ultra performance and ultra value!

### **Convection Cooled**

These models use convection cooling to increase the power range of the UP19K family. The UP19K-30H-H5 takes you to 30 W continuous and 45 W intermittent. Using our oversized custom heatsink the UP19K-50L-H5 pushes that to 50 W continuous and 75 W intermittent.

#### Air and Water Cooled

You want the same small package but for high average power. With a few volts to power the fan the UP19K-110F-H5 handles 110 W of laser power in an amazingly compact package. With water cooling, the UP19K-150W-H5 takes 150 watts. With the DI option it is perfect for clean deionized water cooling systems.

#### Calorimeter Mode

With this option every member of the family can be equipped to measure single shot pulse energies as well as average power. From 20 mJ up to 5 J Q-switched or 25 J long pulse.

## Fiber Optic Option

Optional fiber adapters are available for these detectors.

#### ULTRA SERIES UP12E SPECIFICATIONS

COMMON APPLICATIONS

· Low energy OEM

TYPICAL LASERS

CO.

· High repetition rate YAG (various) 38 -32-012 Ø12 ø12 Photolithography Diode bars Medical Excimer 0 0 0 Ti:sapphire 38 38 Ruby (long pulse) 0 0 0 65 All dimensions in mm 20H MEASUREMENT CAPABILITY Spectral Range 0.19 - 11 µm 0.19 - 11 µm 0.19 - 11 µm Maximum Measurable Power 10 W 20 W 70 W Minimum Detectable Power<sup>a</sup> 1 mW 1 mW 1 mW Rise Time (nominal) b 0.3 sec 0.3 sec 0.3 sec Sensitivity c,d 0.6 mV/W 0.6 mV/W 0.6 mV/W Calibration Uncertainty e ±2.5 % ±2.5 % ±2.5 % Repeatability ±0.5 % ±0.5 % ±0.5 % **Power Resolution** ±0.5 % ±0.5 % ±0.5 % Calorimeter mode Sensitivity 0.95 mV/J 0.95 mV/J 0.95 mV/J Maximum Measurable Energyf 5 J 5.7 5 J Minimum Measurable Energy 0.2 J 0.2 J 0.2 J Minimum Repetition Period 1.5 sec 1.5 sec 1.5 sec Maximum Pulse Width 50 ms 50 ms 50 ms Accuracy with energy calibration option ±5 % ±5 % ±5 % Beam Size Dependence 9 ±0.5 % ±0.5 % ±0.5 % DAMAGE THRESHOLDS Max Average Power (continuous) 10 W 20 W 70 Wh Max Average Power (2 minutes) 15 W 30 W 90 W Maximum Average Power Density<sup>1</sup> 36 kW/cm<sup>2</sup> 36 kW/cm<sup>2</sup> 36 kW/cm<sup>2</sup> Pulsed Laser Damage Thresholds Max Energy Density Peak Power Density 1.064 µm, 360 µs, 5 Hz 5 J/cm<sup>2</sup> 14 kW/cm<sup>2</sup> 1.064 µm, 7 ns, 10 Hz 1 J/cm<sup>2</sup> 143 MW/cm<sup>2</sup> 532 nm, 7 ns, 10 Hz 0.6 J/cm<sup>2</sup> 86 MW/cm<sup>2</sup> 248 nm, 26 ns, 10 Hz 0.3 J/cm<sup>2</sup> 43 MW/cm<sup>2</sup> PHYSICAL CHARACTERISTICS **Effective Aperture Diameter** 12 mm Ø Absorber High Damage Threshold - H5 **Dimensions** 38H x 38W x 14D mm 38H x 38W x 45D mm 38H x 38W x 32D mm Weight (head only) 0.13 kg 0.15 kg 0.19 kg **Effective Area** 1.13 cm<sup>2</sup> 1.13 cm<sup>2</sup> 1.13 cm<sup>2</sup>

105

20H

70W

Beam centered.
Minimum cooling flow 0.5 liter/min, water temperature ≤ 22°C, 1/8 NPT compression fittings for 1/4 inch semi-rigid tube. Contact Gentec-EO for clean deionized water cooling module option.
At 1064 nm.

Specifications subject to change without notice.



GENTEC ELECTRO-OPTICS INC. www.gentec-eo.com



a. Nominal value, actual value depends on electrical noise in the measurement system. b. With Gentec-EO TPM 300CE, DUO, SOLO or P-LINK monitor. c. Maximum output voltage = sensitivity x maximum power. d. Higher sensitivity with internal circuit. Contact Gentec-EO, e. Including linearity with power. With Gentec-EO monitor.

f. For 360 µs pulses. Higher pulse energy possible when customized for long pulses (ms), less for short pulses (ns).





PH100-Si

# PHOTO DETECTORS

**Low Power** 

 $\epsilon$ 

- Fast
- Small Size
- Si, Si-UV, Ge
- Full NIST-Traceability
- Smart Interface



## PH Series Photo Detector Specifications

Gentec-EO offers a series of photo detectors that can extend your laser power measurement capability down to the picowatt range. We have a solution for your low power measurements for wavelengths from 200 nm to 1.8 micron. You have 3 models to choose from: silicon, UV-enhanced silicon and germanium. The silicon detectors provide a generous aperture for laser beams up to 10 mm in diameter. The threaded aperture allows you to install filters, attenuators, fiber optic adapters or other optics that suit your specific needs. All of the detectors have a NIST traceable calibration and Personal wavelength correction™. Our photo detectors are compatible with the Gentec-EO SOLO PE and P-Link monitors.

#### Silicon

This is the choice for typical low power CW lasers from 300 to 1100 nm. Make measurements down to 600 pW and up to 0.75 W with an OD-2 attenuator.

### Silicon-UV

For low power CW lasers working at shorter wavelengths we offer this photo detector sensitive to 200 nm. Maximum power with OD-2 attenuator is 5.4 mW. Other performance characteristics are similar to the regular silicon detector.

#### Germanium

For longer wavelengths, our germanium photo detector is a nice solution between 0.8 and 1.8 microns for laser beams up to 9 mm in diameter. With the attenuator it will go to 0.75 W.

#### Attenuators

To extend the performance of the photo detectors to higher measurable powers you can select one of our attenuators. The transmission of the OD-1 is 10% and the transmission of the OD-2 is 1%.

Diode

Ge

· Raman shifted YAG

Si-UV

HeCd

Wafer inspection • Data storage

Multimode fiber • Bioluminescence


MicroYAG

HeNe

Diode

Argon Ion

73



All dimensions in mm

Control of the Contro	PH100-Si PH100-Si <sup>UV</sup>		PH78-Ge	
FACUREMENT CARARY TO				
EASUREMENT CAPABILITY				
ectral range	300 -1100 nm	200 -1100 nm	800 -1800 nm	
aximum Measurable Power (1064 nm)	30 mW	54 μW	27 mW	
With OD-1 attenuator	0.3 W	0.54 mW	0.27 W	
With OD-2 attenuator	0.75 W	5.4 mW	0.75 W	
inimum Measurable Power <sup>b</sup>	600 pW (980 nm)	600 pW (850 nm)	20 nW (1550 nm)	
With OD-1 attenuator	6 nW	6 nW	200 nW	
With OD-2 attenuator	60 nW	60 nW	2 μW	
se Time (nominal)	0.2 sec	0.2 sec	0.2 sec	
ak Sensitivity	0.5 A/W (980 nm)	0.5 A/W (850 nm)	0.9 A/W (1550 nm)	
libration Uncertainty	±6.5% 0.3-0.44 μm	±6.5% 0.2-0.44 μm		
	±2.5% 0.44-1.0 μm	±2.5% 0.44-1.0 μm	±3% 0.8-1.8 μm	
	±3.5% 1.0-1.1μm	±3.5% 1.0-1.1 μm		
ccuracy (with OD-1 and OD-2 attenuators)	±5%	±5%	±5%	
avelengths available	325, 351, 440, 488, 514.5,	248, 266, 325, 351, 440,	850, 980, 1064, 1310	
rith SOLO PE and P-LINK)	532, 632.8, 670, 694, 780,	488, 514.5, 532, 632.8, 670,	32.8, 670, 1550, 1620 nm	
	850, 980, 1064 nm	694, 780, 850, 980, 1064 nm		
With OD-1 attenuator	440, 488, 514.5, 532,	488, 514.5, 532, 632.8,	980, 1064, 1310,	
	632.8, 670, 694, 780,	670, 694, 780, 850,	1550, 1620 nm	
	850, 980, 1064 nm	980, 1064 nm		
With OD-2 attenuator	632.8, 670, 694, 780, 850,	632.8, 670, 850,		
	980, 1064 nm	980, 1064 nm	N/A	
AMAGE THRESHOLDS				
aximum Average Power Density (1064 nm)	100 W/cm <sup>2</sup>	100 W/cm <sup>2</sup>	100 W/cm <sup>2</sup>	
turation Level (1064 nm)	30 mW/cm <sup>2</sup>	50 μW /cm²	< 100 mW/cm <sup>2</sup>	
HYSICAL CHARACTERISTICS				
fective aperture diameter	11.28 mm Ø	11.28 mm Ø 10 mm		
mensions	36 mm Ø x 34.1 mm	36 mm Ø x 34.1 mm 36 mm Ø x 34.1 mm		
eight	71 g	71 g	71 g	
tive area	1 cm²	1 cm <sup>2</sup>	0.78 cm <sup>2</sup>	

With Gentec-EO SOLO PE and P-Link monitors. Half hour warm-up before offset nulling, offset nulling on the lowest scale after each new power supply. Temperature ±0.5 degrees.

Specifications subject to change without notice.